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*ZLUV856 adopts UV luorescence lamps to simulate sunlight effect. Its inner temperature and humidity can be properly controlled to obtain a continuing condensation on the sample for fully evaluating the damage caused by sunlight, humidity and temperature (effects caused by aging of the material may be fading, intensity reduction, cracking, flaking, chalking and oxidation).*

*The UV light accelerated aging test chamber simulates in days or weeks the damage which, in real conditions, would occur in months or years.*

*UV light can emulate the effect of sunshine, while condensation and water spray system can emulate the effects of rain and dew.*

*It is possible to choose UVA or UVB lamps.*

*UVA lamps provide a true simulation of sunlight in the short-wave range and are the most used to test the effects of sun damage. While UVB lamps, having an even shorter wavelength, simulate a much more severe and less realistic environment and are suitable for testing very durable materials, usually exposed to intense sunlight and for a long period of time.*

*A typical test cycle generally carries out under strong irradiation of UV light or in a dark and wet condensation period with 100% relative humidity.*

*These tests are generally employed in some specific fields, such as paint and coatings, automotive industry, plastic, wood, etc.*

Modello / Model	ZLUV856
Standard	ISO 16474-3 ; ISO 4892-1 ; ISO 4892-3 ; ASTM D 4799 ; ASTM D 4587 ; ASTM D 4329 ASTM G 151 ; ASTM G 154
Lampade UV UV lamps	UVA or UVB 40W x 8 pcs (life about 6.000 hours)
Lunghezza d'onda Wavelength	340 nm (UVA Lamps) or 313 nm (UVB Lamps) *
Campo di misura dell'irradiazione Irradiance range	0.1 W/m <sup>2</sup> ~ 1.55 W/m <sup>2</sup>
Regolazione temperatura del pannello nero Adjustment of black board temperature	+10 ~ +80°C
Area di esposizione campioni Test panel exposure area	5.175 cm <sup>2</sup>
Consumo d'acqua Water consumption	7Lt/day (condensation) 4Lt/minute (spray)
Capacità della camera Chamber capacity	48 pcs (75x150 mm)
Alimentazione / Power supply	220 V
Potenza / Power	3KW
Dimensioni (LxPxA) Size (LxWxH)	1360 x 520 x 1310 mm
Peso / Weight	160 Kg

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## STANDARD SUPPLY

- *Main unit*
- *User's manual*
- *UVA or UVB lamps (to choose)*

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
Variations in light intensity and irradiance can affect the test and degrade the material under test in different ways.

As the lamps wear out they have a normal decay by varying their irradiance. The instrument is equipped with an automatic irradiance control system that compensates for this decay by maintaining a constant irradiation by varying the power of the lamps.

However, prolonged use inevitably leads to continuous decay of the lamps. At this point the operator can calibrate the instrument using a radiometer. If the instrument fails to reach the set point following calibration, the operator must replace the lamps corresponding to the relevant sensors and calibrate again.

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The radiometer is used for the calibration of UVB or UVA lamps. As for UVB lamps, these are correctly calibrated at a wavelength of 313 nm. As for the UVA lamps, they are calibrated at a wavelength of 340nm.

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The spray function reproduces thermal shocks caused by temperature excursions such as rain. Some particular conditions, such as sudden temperature rises due to sunny days, can cause thermal shocks as the temperature of the material undergoes drastic changes. These sudden increases in temperature allow us to evaluate the properties of many types of material. Thanks to the spray function it is possible to reproduce these particular atmospheric conditions.



The instrument is equipped with 12 nozzles, 6 for each side of the test chamber. With the spray function you can set the duration of the spraying times and intervals. These intermittent sprays of water cool the samples rapidly, creating continuous thermal shocks.

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In many outdoor environments, materials are placed in wet condition for over 12 hours each day. Studies have shown that the main factor of this wet condition outdoor is caused by dew, not rain.

The instrument simulates the outdoor moisture erosion through the reproduction of a realistic condition. During the condensation cycle test, water on the bottom of the chamber is heated to obtain superheated steam filling the test chamber. Hot steam makes the chamber maintain 100% relative humidity, and maintain a relatively high temperature. Sample was fixed on the wall of test chamber. One side of the sample surface is exposed to the ambient air of test chamber. While the other side is exposed to the natural environment which has a cooling effect, both internal and external surfaces of the sample have different temperature. This temperature difference leads to the test surfaces always have drips caused by condensation process.

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In each cycle, the temperature can be controlled at a set value. At the same time, the black panel thermometer can monitor the temperature. Temperature increase can accelerate the aging process, and the temperature control is also very important for reproducibility of the test. Black panel thermometer consists of a PT100 sensor and a metal panel painted by black coating, and be exposed to the same condition as test panels. It's used to monitor exposing test panels surface temperature during test. This one can be set at any value according to different requirements and also be controlled automatically during the whole test. It also can be calibrated periodically.

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The test program can be modified independently to conform to the required standards. Each program can be set in 10 segments and each segment can be programmed with 4 different types of work: irradiation, nebulization, condensation, final cycle. It is possible to save up to 6 predefined programs.

The screenshot displays a control interface with a dark blue background. On the left, there are six buttons labeled 'Program 1' through 'Program 6'. The main area shows configuration parameters for a selected program (Program No. 1):

- Segment: 1
- Work Condition: 1
- Time: 4 Hour(Minute)
- Irradiance: 6.77 W/m²
- Temperature: 60 °C
- Set Program Time: 1000 Hour 0 Minute
- Past Program Time: 0 Hour 0 Minute

Buttons for 'Preview', 'Save', 'Clear Time', and 'Return' are visible on the right side of the interface.

The screenshot shows a control interface with a dark blue background. At the top, there are buttons for 'Login', 'Exit', 'Manage Parameter', and 'User: admin'. The main area is divided into several sections:

- Set Parameter:** Delay Hot Fan (5), Door Switch (red), Max. Water Tem. (60.0), Min. Water Tem. (50.0).
- Cumulative Lamp Irradiate Time:** Four rows for Group 1# Lamps, each with a 'Clear' button.
- Manual Control:** Manual Spray, Manual Hot Fan (3200), Manual UV (6.77 W/m²), Manual Water Heat.
- Lamps select:** UVA, UVB (selected), UVC.
- Bottom Status:** Cumulative UV Run Time: 0 Hour 0 Minutes, Manufacturing Code: 111002, and a 'Return' button.