Quick Afla M1 Strip Test - Cod. ASTM1/96 96 Test - Kit immunocromatografico

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Quick Afla M1 Strip Test - Code ASTM1/96 96 Tests - Lateral-flow test kit

Rapid immunochromatographic test for the qualitative/semiquantitative determination of Aflatoxin M1 in milk according to European limits (50 ppt)

Introduction

Quick Afla M1 Strip Test is a rapid assay based on a highly specific immunologic reaction which develops between the anti-Aflatoxin M1 monoclonal antibody of the lateral flow test and the Aflatoxin M1 of the milk sample.

Aflatoxins are toxic and carcinogenic metabolites produced by some moulds such as *Aspergillus flavus* and *Aspergillus parasiticus*. Crops may be contaminated by one or more of the four following sub-types of aflatoxins: B1, B2, G1 and G2. Aflatoxin B1 is the most toxic and frequently detected type.

When cows eat contaminated feed, Aflatoxin B1 is converted by hydroxylation to Aflatoxin M1, which is subsequently secreted in the milk and urine of lactating cows. It is the most hepatic carcinogenic agent for some animal species. Since young animals are more sensitive to the toxic action of all aflatoxins when compared to adults, great attention is paid to the analysis of Aflatoxin M1 in milk, as it is widely consumed by babies. As a matter of fact, Aflatoxin M1 is quite stable towards the normal milk processing methods such as pasteurization and, if present in raw milk, it may persist into final products for human consumption. Furthermore, the use of contaminated milk for the production of dairy products does not reduce or destroy the toxin, therefore the determination of this substance at every level of milk distribution is very important.

Most controlling government agencies worldwide have regulations regarding the amount of aflatoxins allowable in human and animal foodstuffs. Many countries established limits for the presence of Aflatoxin M1 in milk and milk products. Since 1998, in the EU the limit for the presence of M1 in milk and reconstituted milk powders has been set at 0.05 mg/L or 50 parts per trillion (50 ppt). CODEX Alimentarius and FDA have established a maximum level of 0.5 mg/L or 500 parts per trillion (500 ppt).

The use of this kit is meant for the determination of Aflatoxin M1 in milk samples in a few minutes and without any need of incubators or additional chemicals/consumables. The presence of preservatives in the samples does not alter the test results.

Principle of the test

Quick Afla M1 Strip Test is a competitive immunochromatographic assay based on the lateral flow format. Each test strip includes a highly specific anti-Aflatoxin M1 monoclonal antibody linked to colloidal gold particles. The milk sample (200 µL) is added to a clear plastic reaction microwell, and used to resuspend the lyophilized reagents to a uniform pink color in the bottom of the microwell. The milk is incubated briefly (5 min) to allow the anti-Aflatoxin M1 antibody bound to the gold particles to engage with the Aflatoxin M1 present in the milk sample. The test strip is then put into the sample well with the arrows pointing downward so that the milk is capillary absorbed by the strip. Any antibody linked to gold particles that is not complexed with the Aflatoxin M1 present in the milk will bind to the Aflatoxin M1 imprinted at the Test line (T-Line), forming a signal (red line) at that position. If the antibody linked to the gold particles has engaged with the Aflatoxin M1 present in the milk sample, the gold particles will flow past the T-Line and reach the Control line (C-Line).

Visual interpretation of the test results: if the T-Line signal intensity is stronger than the signal at the C-Line it means that the result is lower than 50 ppt. If the signal at the T-Line is less intense compared to the C-Line it means that the concentration of Aflatoxin M1 in the milk sample is higher than 50 ppt. If the two lines have equal signal intensities and it is not possible to distinguish which one is more intense by eye, then the test result is about 50 ppt. When only the C-Line is visible and the T-Line cannot be seen, the test result is strongly positive (*equal or higher than 300-400 ppt*).

The kit includes all that is needed to perform the test on milk samples everywhere, on every kind of table or bench.

Precautions

1. Store all components of the kit at 2°-8°C and far from the light. DO NOT FREEZE. Do not use the kit after its expiry date printed on the external label. The kit shelf life is 12 months after the production date, when correctly stored.
2. Before opening the test strips container, allow it to reach the room temperature (25°C ± 5°C).
3. By means of scissors, cut and remove the exact number of reaction wells the remove the needed number of test strips according to the number of samples to analyze. Close back all the unused wells and strips in their container, then promptly put them back in the refrigerator.
4. Follow the instructions as given in this document, strictly. Any variation to times and procedures may cause wrong results.
5. Follow the Good Laboratory Practice when the kit is used and during the analysis. Handle the test strips with care with clean and dry hands.
6. It is suggested to include a negative control sample in every test session. This is not mandatory, but it may be useful to check the kit is working properly.

7. Mix well all milk samples before testing them. Preferably, analyze all samples at room temperature ($25^{\circ}\text{C} \pm 5^{\circ}\text{C}$), however the temperature of the samples does not affect the test results, unless they are still partially frozen.
8. Use a new plastic tip for each sample. Obviously, the reaction microwells and the test strips cannot be reused.
9. To check and follow the reaction times correctly, it is strongly suggested not to run more than 8 tests simultaneously. Start the incubation time at the beginning of each period and for every sample, appropriately, to avoid different incubation times from sample to sample.
10. Quick Afla M1 Strip Test is a screening test: all results should be confirmed by the reference method (HPLC).
11. As suggested by Good Laboratory Practice, it is strongly recommended to wear protective gloves and all usual safety devices when the kit is used and during the analysis.
12. Waste all materials, containers, reagents and samples (and every product that gets in touch with them) appropriately after their use. Do not pour the unused reagents back in their original containers.
13. The presence of preservatives in the samples, such as Bronopol or Azidol, sodium azide, does not affect the test results. It means that preserved samples can be analyzed without any problem.

Test procedure

Preparation of the samples

Liquid milk from every process (fresh, pasteurized, UHT, microfiltered, raw, whole, partially skimmed, skimmed, etc.):

Use the sample as it is, without any pre-treatment. It is not necessary to defat it before the test.


Powdered milk:

Dissolve the sample in distilled or deionized water according to the appropriate proportion weight/volume, as described by the producer of the powdered milk. Shake well the dissolved sample in order to be sure about its perfect solubilization and the absence of residuals of powder in the solution.

Use such a solution as a liquid milk sample (*see the previous paragraph*).

Analysis


1. Take the kit out of the refrigerator.
 2. Before opening the container of the test strips, allow it to reach room temperature ($25^{\circ}\text{C} \pm 5^{\circ}\text{C}$).
 3. Using scissors, carefully cut and remove the right number of wells and test strips according to the number of samples that will be tested. Close back all the unused wells and strips in their container, then promptly put them back in the refrigerator. Open the reaction microwells and insert them firmly in the supplied white plastic frame.
 4. To check and follow the reaction times correctly, it is strongly suggested not to run more than 8 tests simultaneously. Start the incubation time at the beginning of each period and for every sample, appropriately, to avoid different incubation times from sample to sample.
 5. Mix well all milk samples before testing them. Preferably, analyze all samples at room temperature ($25^{\circ}\text{C} \pm 5^{\circ}\text{C}$), however the temperature of the samples does not affect the test results, unless they are still partially frozen.
 6. Aspirate 200 μL of milk by using the supplied 200- μL micropipette with a new plastic tip firmly inserted. Be sure to take this volume when the tip enters deeper than 1 centimeter below the higher level of the milk sample.
 7. Dispense this volume in a new reaction microwell, fixed in the provided frame, being sure to depress the plunger to completely expel the milk sample from the plastic tip into the well.
 8. Using the same pipet tip, **aspirate the sample up and down about 8-10 times to completely resuspend the lyophilized gold particles in the milk sample**. Avoid the production of bubbles. The sample should turn a uniform pink color. After resuspending the particles, remove and discard the pipet tip.
- Note:** if this step is missed or not performed correctly, the test results may be affected.
9. **Incubate the sample for 5 minutes at room temperature** ($25^{\circ}\text{C} \pm 5^{\circ}\text{C}$).
 10. With clean and dry hands, insert the bottom of the test strip into the microwell containing the milk sample. The test strip should be inserted with the printed arrows pointing down. Be sure the test strip is oriented vertically and it reaches the bottom of the well.
 11. **Incubate the test strip for 5 minutes at room temperature** ($25^{\circ}\text{C} \pm 5^{\circ}\text{C}$).
 12. Take the test strip out of the microwell, put it on a horizontal surface and **interpret the result within 5 minutes**, by eye or by means of the optional RDS-1500 PRO Strip Reader (*see the following paragraph*). When using the RDS-1500 PRO Strip Reader for the quantitative interpretation, in order to grant the maximum accuracy and reproducibility of the results, please read the strip test immediately after the 10-minute incubation period (5' + 5') is over.

Interpretation of the results

At the end of the last incubation time, observe the test strip and the red lines that have appeared on it. The line located close to the arrows (that is, close to the zone that has been inserted in the microwell) is called T-Line (Test Line); the line located far from the arrows is called C-Line (Control Line).

- **The Test Line (T-Line) looks more intense than Control Line (C-Line):** the concentration of Aflatoxin M1 in the analyzed sample is lower than 50 ppt.
- **The Test Line (T-Line) shows the same intensity as the Control Line (C-Line):** the concentration of Aflatoxin M1 in the analyzed sample is equal to approx. 50 ppt.
- **The Test Line (T-Line) looks less intense than Control Line (C-Line):** the concentration of Aflatoxin M1 in the analyzed sample is higher than 50 ppt.
- **The Test Line (T-Line) is absent and only the Control Line (C-Line) is visible:** the concentration of Aflatoxin M1 in the analyzed sample is very high (equal to or higher than 300-400 ppt).
- **If no line is present in the C-Line area:** the test is INVALID. It is suggested to repeat the test from the start.

Note: Quick Afla M1 Strip Test also gives an approximate indication of the concentration of Aflatoxin M1 in the sample according to the intensity level of the two lines. As an example, when the Test Line (T-Line) is more intense than Control Line (C-Line) and this one is scarcely visible, the weaker the C-Line looks and the closer to zero the concentration of Aflatoxin M1 in the sample is; similarly, when the Test Line (T-Line) is less intense of the Control Line (C-Line) and the T-Line is scarcely visible, the weaker the T-Line looks and the higher than 50 ppt the concentration of Aflatoxin M1 in the sample is. Such a semiquantitative interpretation is always approximate if realized by eye: for this reason, it should be confirmed by a quantitative method (ELISA or HPLC).

It is possible to keep the test strip and its test result for future references: to this purpose, cut the lower part of the strip (where the arrows are printed) by means of scissors to avoid that the sample flow keeps on running along the test strip and the test result changes.

The lack of liquid migration along the test strip brings to an invalid test. This problem may occur when the test is performed with abnormal milk, coagulated milk or when solid residues are present, or even if the test procedure has not been followed correctly. If such an inconvenient occurs, please repeat the test from the start. All samples that have been interpreted as positive should be confirmed by the reference method (HPLC).

The test strip can be also objectively read by means of the RDS-1500 PRO Reader, which provides quantitative results (in ppt). To start the reading, insert the strip in the black cartridge adapter in the direction as indicated in the manual of the reader. Then, by handling the semi-circular part of the adapter with the strip facing the equipment, put the adapter into the reading chamber located on the upper part of the equipment. Press the "Run test" key on the touch-screen: the quantitative result will be displayed immediately. Please, read the user manual of the RDS-1500 PRO Reader for more information about all other functions and operations featured by the equipment. **When using the RDS-1500 PRO Strip Reader for the quantitative interpretation, in order to grant the maximum accuracy and reproducibility of the results, please read the strip test immediately after the 10-minute incubation period (5' + 5') is over.**

Kit content

- 96 reaction microwells and 96 test strips, parted in cylindrical tubes with drying caps. Each cylindrical tube includes 8 microwells and 8 test strips.
- 200- μ L micropipette, reusable
- 96 disposable plastic tips
- 96-hole reusable frame, to keep the reaction microwells in the appropriate vertical position
- Test manual

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