



UltraSnap™ | Surface ATP Test

UltraSnap is a user-friendly ATP surface test designed for use with Hygiena luminometers. This all-in-one test contains a pre-moistened swab bud for optimal sample recovery. A unique liquid-stable reagent that provides superior reproducibility and robustness and Hygiena's patented Snap-Valve™ technology for easy activation.

Benefits:

- Pre-wetted swab maximizes sample recovery and breaks through biofilm
- Unique liquid-stable reagent gives superior sensitivity and reliable results
- 100% recyclable plastic
- 15 month shelf life at refrigerated temperatures (2 -8 C)
- 4 week shelf life at room temperatures (21 -25 C)
- Costs up to 50% less than other ATP test devices



Watch an instructional demo at youtube.com/HygienaTV



Find support documents, instructional videos, and more at www.hygiena.com


Sensitivity:

Hygiena System	UltraSnap Detection Limit				
femtomoles ATP	10 fmol	1fmol	.1 fmol	.01 fmol	.001 fmol
ATP (mol)	1×10^{-14}	1×10^{-15}	1×10^{-16}	1×10^{-17}	1×10^{-18}
SystemSURE Plus	●●●●				
EnSURE	●●●●●●				
Pi-102	●●●●●●●●●●				

Performance:


UltraSnap's superior performance was validated in a third party study conducted by a large reference laboratory. In addition to superior repeatability (Diagram 1), Hygiena's system demonstrated outstanding performance for linearity, sensitivity, and accuracy. Learn more at www.hygiena.com or contact Hygiena for a full copy of the study.

Repeatability: Variability in results with the same sample



Catalog No.	Description	Quantity
US2020	UltraSnap Surface ATP Test	100

UltraSnap Procedure:



Related Products:



SystemSURE Plus

SystemSURE Plus is used for ATP hygiene monitoring. This palm-sized instrument is easy to use, extremely sensitive and affordable, allowing users to quickly determine the cleaning efficiency and hygienic status of surfaces and liquid samples to ensure product quality and reduce costs.



EnSURE

EnSURE is a high-sensitivity luminometer that measures multiple quality tests such as ATP, TVC, Coliform, *E. coli* and more all on one system. The monitoring system includes SureTrend data tracking software.



AquaSnap

AquaSnap measures ATP found in liquid samples, utilizing a honeycomb dipper to consistently collect 100 µl of liquid sample. The device uses a liquid-stable reagent in place of freeze-dried enzymes giving it superior reproducibility and sensitivity. AquaSnap is available in Total and Free format.



SuperSnap

SuperSnap is the most sensitive ATP surface test in Hygiena's ATP hygiene monitoring portfolio of all-in-one devices. Designed to work with all Hygiena luminometers, SuperSnap is an all-in-one test which detects extremely low levels of ATP. SuperSnap is predominantly used in food and beverage processing facilities as an allergen cross contamination prevention.

Catalog No.	Description	Quantity
SS3	SystemSure Plus ATP Monitoring System with SureTrend Software	1
ENSURE	EnSURE Monitoring System with SureTrend Software	1
AQ-100X	AquaSnap Free Water ATP Test	100
AQ-100FX	AquaSnap Total Water ATP test	100
SUS3000	SuperSnap	100