

Evaporatore Rotante LabTech Ev

Modelli Disponibili

Technical Specifications

Model	EV311-V	EV311-V Plus	EV311-L	EV311-L Plus
Rotation Speed Range	20 to 340 rpm	20 to 340 rpm	20 to 340 rpm	20 to 340 rpm
Rotation Speed Display	LED	LED	LED	LED
Adjustable Head Angle	0° to 60°	0° to 60°	0° to 60°	0° to 60°
Motor Power	40 Watts	40 Watts	40 Watts	40 Watts
Stroke Displacement	160mm Motorized	160mm Motorized	160mm Motorized	160mm Motorized
Unit Safety Protection	Automatic lift up when electric power shut down	Automatic lift up when electric power shut down	Automatic lift up when electric power shut down	Automatic lift up when electric power shut down
Temperature Range	Ambient to 180°C	Ambient to 180°C	Ambient to 180°C	Ambient to 180°C
Material	304 SSL	PTFE Coating	304 SSL	PTFE Coating
Temperature Controller Display	LED	LCD	LED	LCD
Temperature Controller Mode	P.I.D.	P.I.D.	P.I.D.	P.I.D.
Temperature Accuracy	±1°C (Water) ±2°C (Oil)	±1°C (Water) ±2°C (Oil)	±1°C (Water) ±2°C (Oil)	±1°C (Water) ±2°C (Oil)
Heating Power	1300 Watts	1300 Watts	1300 Watts	1300 Watts
Bath Diameter	250 mm	250 mm	250 mm	250 mm
Bath Capacity	5 Liters	5 Liters	5 Liters	5 Liters
Evaporation Volume Range	500 to 3000 ml	500 to 3000 ml	500 to 3000 ml	500 to 3000 ml
Bath Safety	Automatic stop	Automatic stop	Automatic stop	Automatic stop

Distribuito da:
Zetalab s.r.l.



Via Castelfidardo, 11 - 35141 Padova
Telefono 049 2021144 - Fax 049 2021143
Internet: www.zetalab.it - e-mail: info@zetalab.it