

Telecamera Wi-Fi: **WF5**



Wi-Fi Camera

WF5

Nuova telecamera WiFi, con attacco a passo C.

Sensore da 5Mp per condividere le immagini in tempo reale al meglio.

Trasmissione WiFi standard, semplice da trovare e da utilizzare con qualsiasi dispositivo (cellulare-tablet ecc)



Caratteristiche tecniche:

Risoluzione:	5 MP
Sensore:	1/2.5" CMOS True color
Risoluzione:	2592 X 1944
Range di frequenza:	2.4000 – 2.4835 GHZ
Modalità di utilizzo:	WiFi o USB
Esposizione:	Manuale o automatica
Bilanciamento del bianco:	Manuale o automatica
Compatibilità browser:	Safari – Chrome – I.E.
Frame rate:	JPG Output: 640x480@40fps,800x600@40fps,1280x960@40fps,1280x720@40fps,2592x1944@15fps
Transmit:	Bandwidth 150M Bps
WiFi:	IEEE 802.11b/g/n
Alimentazione:	Separable and rechargeable lithium battery (3.7V 1900mAh) ; DC 5V USB alimentatore per PC ; External DC 5V 1.5A adapter
Sistemi operativi:	Microsoft® Windows® XP / Vista / 7 / 8 (32 & 64 bit) MAC OS X (Wi-Fi Working Mode) IOS5/IOS6 (Wi-Fi Working Mode) Andriod 2.2 or above
Requisiti minimi:	CPU: Intel Core2 1.6GHz or Above Memory:2GB or Above USB Port : 2.0 DVD-ROM(Required When Installing Software)
Optical port:	C-Mount

