



Sension+ 5203 pH liquid filled combination electrode with encapsulated Ag/AgCl reference element, screw cap S7.

Codice articolo: LZW5203.97.0002

pH liquid electrode (refillable) with 2 ceramic diaphragm. Glass body. Range 0-14 pH, 0 to 80 °C. S7 screw cap. This pH electrode is suitable for aqueous media in general up to pH 14.

Fast responding.

For samples containing sulphurs, reducing sugars or other substances which react with silver ions.

Ag⁺ ion barrier

Specifiche tecniche

Accuratezza:	0.02 pH with pH3 Meter
Connettore:	Attacco a vite S7
Diametro:	12 mm
Dimensioni (DXL):	12 mm x 130 mm
Giunzione sonda:	2 x Ceramic pin
Intervallo di misura:	0 - 14 pH
Intervallo di temperatura:	0 - 80 °C
lunghezza:	130 mm
Materiale:	Sensor Body: Glass
Parametro:	pH
soluzione di riempimento:	LZW9500.99 KCl 3M
Tipo:	Lab
tipo di sensore:	Glass
Tipo di sonda:	Laboratory