



Rotary Evaporator Line

- High performance and low noise
- Modular concept
- Motorized lift
- LED/LCD display
- PTFE+FKM anticorrosion sealing system
- On-line sample addition
- Support both water bath and oil bath
- Self-calibration function
- Self-protection



Best price/quality ratio in the market

Technical Specifications

Model	EV311-V/L	EV311PLUS-V/L
Main Unit		
Rotation speed	20-240 rpm	
Rotatio speed display	LED	LCD
Head titled angle	0-60°	
Motor power	40 W	
Lifting mode (Motor/Manual)	Motor	
Adjustable height	160 mm	
Self protection	Automatically lift up when no power supply	
Dimensions	380 mmL x 340 mmW x 692 mmH	
IP class	IP20	
Working temperature	0-40 °C	
Weight	23 Kg	

Heating Bath		
Temp. range of heating bath	Ambient temperature to 180 °C	
Heating bath material	204 SSL	Teflon Coating
Temperature control mode	PID - LED Display	PID - LCD Display
Temperature accuracy	±1 °C (Water Bath) / ±2 °C (Oil Bath)	
Heating power	1300 W	
Heating bath diameter	250 mm	
Volume	5 L	
Self-protection	Automatic stop heating when dry bath condenser column	

Condenser Column	
Evaporation volume	50-3000 ml
Condenser type	Vertical or Diagonal
Condenser area	1200 cm ²

Distribuito da:
Zetalab s.r.l.



Via Castelfidardo, 11 - 35141 Padova
Telefono 049 2021144 - Fax 049 2021143
Internet: www.zetalab.it - e-mail: info@zetalab.it

