

LCv 4010

MediLine



Total gross capacity	361 litre
Ext. dimensions in mm (w/d/h)	600/615/2000
Int. dimensions in mm (w/d/h)	440/435/1105
Energy consumption in 24 h	1,800 kWh
Climate rating ¹	SN-ST (+38°C)
Cooling system	static
Defrosting	manual
Temperature range - fridge	+3°C to +16°C
Temperature range - freezer	-9°C to -30°C
Evaporator cover - fridge	yes
Housing material/colour	steel/white
Door/lid material	Steel
Liner Material	moulded polystyrol white
Type of control	electronic controller
Temperature display	external digital
Interior light	–
Fault: warning signal	visual and audible
Interface	RS 485
Volt-free alarm contact	yes
Shelves total	5
Adjustable Shelving	4
Raised base shelf	yes
Shelf material	glass
Shelf loading	24 kg
Door handle	Ergonomic door handle
Self-closing door	yes
Insulation	70 - 70 mm
Refrigerant	R 600a
Noise level	52 dB(A)
Gross/net weight	91/85 kg
EAN-Nr.	9005382181751
Connection Rating	1.5 A / 200 W

¹ Climate Ratings: **SN** +10°C to +32°C **N** +16°C to +32°C **ST** +16°C to +38°C **T** +16°C to +43°C