



pH electrode, HF resistant (<1 g/L)

Codice articolo: CRI5206.99

Combined pH electrode with a specific membrane for use in samples containing Fluoridric Acid with a HF concentration <1 g/L and a pH ≥ 3 . Cylindrical membrane, diaphragm reference. Not for use in viscous solutions, for samples with colloidal particles, or distilled water.

Distribuito da:
Zetalab s.r.l.

Via Umberto Giordano, 5 - 35132 Padova
Tel 049 2021144 - Fax 049 2021143
www.zetalab.it - email: info@zetalab.it

Zetalab.it