



## Sension+ pH combination electrode (dirty samples)

Codice articolo: LZW5011.97.0002

pH combination electrode with open circular diaphragm and solid electrolyte. Glass body. Range 2-14 pH, 0-80 °C. 1 m cable. For benchtop meter.

The electrode is suitable for waste water, dirty samples, viscous media, and food applications. The open diaphragm around the membrane facilitates the measuring in some media where other electrodes have certain difficulties.

---

### Specifiche tecniche

Accuratezza:	0.02 pH with pH3 Meter
caratteristica speciale:	Solid gel for high solids content
Connettore:	BNC
Diametro:	12 mm
Dimensioni (DXL):	12 mm x 130 mm
Giunzione sonda:	Open
Intervallo di misura:	2 - 14 pH
Intervallo di temperatura:	0 - 80 °C
lunghezza:	130 mm
Lunghezza del cavo:	1 m
Materiale:	Sensor Body: Glass
Parametro:	pH
soluzione di riempimento:	Non-refillable solid polymer
Tipo:	Lab
tipo di sensore:	Glass
Tipo di sonda:	Laboratory