



Sension+ 5202 pH liquid filled combination electrode with screw cap S7 for general use. Large, semi spheric membrane, very robust and fast responding.

Codice articolo: LZW5202.97.0002

pH liquid electrode (refillable) with ceramic diaphragm. Glass body. Range 0-12 pH, 0 to 80 °C. S7 screw cap. This pH electrode is suitable for aqueous media in general. Fast responding due to large, semi spheric pH membrane.

Fast responding due to large, semi spheric pH membrane.

Specifiche tecniche

Accuratezza:	0.02 pH with pH3 Meter
Connettore:	Attacco a vite S7
Diametro:	12 mm
Dimensioni (DXL):	12 mm x 130 mm
Giunzione sonda:	Ceramic
Intervallo di misura:	0 - 12 pH
Intervallo di temperatura:	0 - 80 °C
lunghezza:	130 mm
Materiale:	Sensor Body: Glass
Parametro:	pH
soluzione di riempimento:	LZW9501.99 KCl 3M + AgCl
Tipo:	Lab
tipo di sensore:	Large, semi-spheric glass membrane.
Tipo di sonda:	Laboratory