

VIBROMIX 50

Large shaker


Introducing our large benchtop shaker, built for continuous operation and heavy workloads. This compact yet robust unit offers a 25 mm orbital shaking motion and features a variety of platforms and clamps to suit your needs. Elevate your lab's shaking capabilities with this reliable and versatile tool.

TECHNICAL SPECIFICATIONS

Model	Vibromix 50
Speed set range	50 rpm to 500 rpm, step 1 rpm
Control	digital
Motion	orbital
Orbit	25 mm
Max. load capacity	10 liters
Available platforms	610 x 460 mm 760 x 460 mm
Display	LED
Timer	0 to 99 min and continuous, steps in 1 sec
Max. power consumption	250 W
Dimension (H x W x D)	150 x 600 x 550 mm
Crating (H x W x D)	225 x 750 x 750 mm
Weight (Net/Crated)	48,0 / 71,0 kg
Motor type	DC brushless motor
Permitted working ambient	"temperature 4°C to 75°C, rh 75% non-condensing"

ORDER INFORMATION

230 V, 50/60 Hz	Vibromix 50
Code	621.000
Order ref.	283271
120 V, 50/60 Hz	Vibromix 50 120 V
Code	1.621.000
Order ref.	392720



For more information visit Domel website: <https://labs.domel.com> / Technical specifications subject to change.

Via Umberto Giordano, 5 - 35132 Padova


Tel 049 2021144 - Fax 049 2021143

www.zetalab.it - info@zetalab.it

Zetalab.it

Distribuito da: **Zetalab s.r.l.**

Platforms and clamps



Sticky plate 610 - platform only for sticky rubber


Code 621.222

Order ref. 301449

Sticky rubber, transparent blue

Code 740.240


Order ref. 462093



Bar plate 610 - Universal platform with squeeze bars

Code 621.300

Order ref. 300186




Uni plate 610 - Universal platform (for clamps)

Code 621.250

Order ref. 283432

Clamps

Nominal volume	125 ml	250 ml	500 ml	1000 ml	2000 ml
No. of clamps/ plate	48	35	20	12	8
Order ref.	169129	169130	204764	204832	204853



Uni plate 760 - Universal platform (for clamps)

Code 621.450

Order ref. 328927

Clamps

Nominal volume	125 ml	250 ml	500 ml	1000 ml	2000 ml
No. of clamps/ plate	60	40	24	15	12
Order ref.	169129	169130	204764	204832	204853



Clamps

Nominal volume	25 ml	50 ml	125 ml	250 ml	500 ml	1000 ml	2000 ml
Order ref.	185520	169143	169129	169130	204764	204832	204853